

LOCAL VIBRATION THERAPY FOR
VETERINARY TREATMENTS





1

PROVEN EFFECTIVENESS

Pain reduction and improved motion thanks to NOVAFON

NOVAFON's gentle vibrations help to reduce pain and alleviate symptoms of neurological disorders.



Increase quality of life



Maintain Thoroughness and mobility



Reduction of tendon pain



Improve gait sequences



Alleviation of muscle pain



Relief of osteoarthritic pain



Alleviation of joint pain



Support mental and physical relaxation



Improve bonds



Preparation and regeneration of muscles (in sports)

2

THE ADVANTAGES OF THE NOVAFON TREATMENT

Gentle yet deep acting — from head to toe for horses and small animals



VERTICAL STIMULATION

Due to its vertical stimulation, the effect unfolds deep in the tissue. In this way, the vibrations also reach deep tissue treatment areas.



GENTLE TREATMENT

NOVAFON's gentle technology enables a wide range of treatments from head to toe. This means that even sensitive areas of the body can be treated in a pleasant manner.



6 CM+ (2.5 INCHES) PENETRATION DEPTH NOVAFON

devices have a penetration depth of 6 cm+, in contrast to other therapy devices that only offer surface stimulation.



Discover our smartest new addition to the NOVAFON family: DAS NOVAFON power 2. You can now enjoy an unequalled range of intensities, as well as 3 frequency levels that make your NOVAFON application even more unique.

3

SAY HI! TO THE FUTURE OF RECOVERY DAS NOVAFON power 2



4


DISCOVER THE NOVAFON APP.

Treat it right:
The NOVAFON application guide

THE PERFECT ADDITION TO YOUR NOVAFON SOUND WAVE DEVICE

The NOVAFON application guide offers you the opportunity to get a better understanding of your device, and to find, understand and carry out the right treatment for you and animals.

Additionally, the app supports you in integrating the NOVAFON better into your everyday life.



- Receive useful information about the device and how it works
- Find out what you need to know about areas of the body and backgrounds
- Provides step-by-step instructions on how to use the device


Scan the QR code & download the app for free



5

THE FUTURE IS SMART Useful Bluetooth connectivity

The direct Bluetooth connection between DAS NOVAFON power 2 and the free NOVAFON app provide you with an even greater and customizable range of settings. Immerse yourself in the smart future of local vibration therapy.



ADJUSTABLE INTENSITY LEVELS

The large range of intensities with 3 frequencies and 15 intensity levels offers the most diverse range of settings to date, individually adaptable for all animal patients.

CUSTOMISABLE TREATMENT

Settings for cats, dogs or horses at the touch of a button: For the most effective use, you can configure individual treatment programmes. You can save these for repeated use in the practice and for your own application. This makes DAS NOVAFON power 2 an all-rounder for all uses.

GET CREATIVE

Make DAS NOVAFON power 2 your very own everyday companion and let it greet you or your pet. Select your favourite display colour.



6

PROVEN TECHNOLOGY

Successfully used for many years by animal therapists, pet owners and professional athletes.

"In every horse that I treated with NOVAFON during the test phase, I noticed a striking relaxation of the treated muscles and a significant warming of the treatment areas."

Animal therapist K.^{1,4}



"In particular, dogs with weak muscles in the hindquarters show significant tension in the lower lumbar area. Conversely, permanent pain in this area can also lead to disturbed movement patterns and even to injuries. So the NOVAFON has now found a permanent place in our preparation and warm-up before exercise."

Veterinary practitioner M. for physiotherapy and naturopathy^{2,4}

"We're really enthusiastic about NOVAFON. I use it for our ponies, for the dog and also for us. We love it dearly! The customer service and the advice is also great and there is always somebody available to help with our problems. Even as a layperson, it was explained to me in an understandable way how to use the NOVAFON on my individual patients."

User B.^{3,4}

¹ Website, Testimonial from Hands 4 Horses: novafon.com/media/pdf/8f/ee/a8/Erfahrungsbericht_Pferd_Hufe_Muskelrelaxion_Hands4Horses.pdf. ² Instagram, Therapist report from TierheilPraxis. ³ Google review from a user ⁴ Case report: There may not be sufficient clinical data to verify the general adequacy of local vibration therapy for the symptoms described. These are some quotes from users:



"It's possible to use the NOVAFON in a very targeted and individual way. This benefits us as we can treat every horse in the stable according to its individual needs. SAP Hale Bob OLD especially enjoys regular treatment for the back and neck."

Ingrid Klimke

Olympic champion and NOVAFON brand ambassador
Ingrid Klimke



Discover many more
testimonials
on Instagram
@novafon_vet

LOCAL VIBRATION THERAPY

Proven effectiveness, variety of applications

THREE FREQUENCY SETTINGS

With a total of three different frequencies, 100, 75 and 50 Hz, the NOVAFON enables each individual area of the body to be treated systematically.

100 Hz

RELIEVING MUSCLE PAIN When the sensory stimulus of vibration interacts with nerve fibres, the transmission of pain signals to the brain can be interrupted as the sensory stimulus overrides the pain stimulus. This effect, known as "Gate Control Theory", can last for a long time after treatment if used regularly.

In this way, pain is alleviated — without medication.

50 Hz

RELEASE OF TENSION

The deep-acting vibrations release muscle tension and loosen hardened tissue from head to foot: From the jaw muscles to the hindquarters.


75 Hz - ONLY DAS NOVAFON POWER 2

MUSCLE ACTIVATION

The gentle vibrations activate the muscles.

ACTIVE REGENERATION

Deep tissue stimulation fosters the body's regeneration mechanisms. This allows the muscles to recover more quickly and enhances physical performance. Blood flow is stimulated, supplying the muscles with more metabolites, such as oxygen.



COMPARISON







Which NOVAFON is best for me?

You can use our comparison section to find out which NOVAFON model is best suited to you and your needs.

With DAS NOVAFON power 2 and its choice of 3 frequencies - each with 15 intensity setting - you can target all areas of the body and treatments.

If you predominantly opt for more gentle treatment for sensitive areas of the body, we recommend DAS NOVAFON.

For treatment of large and compact muscle groups, e.g. for athletes, DAS NOVAFON power, when used at a higher intensity on levels 2 and 3, may be the best choice.

	DAS NOVAFON POWER 2	DAS NOVAFON POWER	DAS NOVAFON
			
USER GROUPS	For personal and professional use. For people & animals – all areas of the body	For personal and professional use. For people & animals – all areas of the body	For personal and professional use. For people & animals – especially sensitive areas of the body
INTENSITY LEVELS			
AREAS OF APPLICATION	Pain, tension, sports injuries, scars, CMD, joint pain/arthrosis, logopaedic & neurological symptoms and regeneration	Pain, tension, sports injuries, scars, logopaedic & neurological symptoms and regeneration	Pain, mild tension, sports injuries, scars, CMD, joint pain/arthrosis, logopaedic & neurological symptoms
TECHNICAL SPECIFICATIONS*			
Frequencies	100 Hz 75 Hz 50 Hz	100 Hz 50 Hz	100 Hz 50 Hz
Amplitude	0.3 – 4 mm	0.5 – 4 mm	0.3 – 3.5 mm
Water resistance	IP44 - Splash-proof	IP22 - Drip-proof	IP22 - Drip-proof
Dimensions	83 x 187 x 38 mm (W x H x D)	83 x 187 x 38 mm (W x H x D)	83 x 187 x 38 mm (W x H x D)
App connectivity	2.4 GHz	-	-
App	Settings & display customisation	-	-
	course of treatment tracking	Treatment videos possible on the app	Treatment videos possible on the app
	default settings optimisation	-	-


*detailed technical information can be found in the product data sheets in the download area of our website <https://novafon.com/en/downloads>

14

EASILY INTERCHANGEABLE ATTACHMENTS

NOVAFON has an intelligent system for connecting attachments.

All attachments can be easily connected/disconnected by rotating them one quarter turn.



Disc attachment
Treatment of larger body regions in a sweeping motion



Ball attachment
For spot treatment of smaller areas

Disc and ball attachment, as well as charging cable and power plug, are included when you buy a NOVAFON device.


<https://novafon.com/en/products>

15

ADDITIONAL ATTACHMENTS

Optimise your treatments with these attachments*


Rehabilitation set
Disk attachment XL ideal for treating larger muscle groups
Pin attachment particularly Stimulating Cone attachment for precise spot application



Magnetic attachment
Deeper spot treatment of, for example, trigger points or areas located deep beneath the skin



Sensory roller set
Combines the blood flow stimulation of the spiky ball massage with the vibrations of the NOVAFON sound wave device, and stimulates surface sensibility as well as deep sensibility



Physio set
Fascia ball for treatment on concave body areas
Fascia fork enables the application of spot pressure, for example next to the spine
Fascia wedge for sweeping large muscle areas and deep-set fascial tissue



*All attachments are available in black, anthracite gray and white.

16

MODULAR SET

The first set of attachments developed particularly for use on animals

Trigger module effective bilateral stimulation, e.g. intercostal muscles.

Grip module suitable for fascia treatment and places sensory input on the skin.



Oval attachment is ideal for wider treatment on animals.

Nub module is particularly stimulating and glides smoothly through fur.

17



The new 4-in-1 modular set with the different surface profiles unites two effective treatment principles:

Blood flow stimulation with the attachment surfaces as well as **pain alleviation and muscle stimulation** from the NOVAFON vibrations. The smart click system makes it easy to change the attachment during treatment. After use, clean the attachment with soap under running water.

The modular set comes in an exclusive and practical transport case with belt loop – perfect for mobile applications, in stables or on the go.



18 ACCESSORIES

For the best NOVAFON experience



With the practical extension handle you can simply work on large animals without the need for a stool. First, insert the NOVAFON device into the extension handle, lock into place and continue treatment as normal.



In order to benefit from your NOVAFON products for a long period of time, we offer you the best way to clean them in a gentle but hygienic way: NOVACLEANER. Completely alcohol-free, developed for medical products, dermatologically tested and particularly gentle on the skin. Available individually or as part of a set with our soft microfiber cloth.



All attachments are available in black, anthracite grey and white.

19 PRACTICAL CHARGING STAND

Simplified charging and professional presentation



The charging stand complements the NOVAFON device perfectly: Its outstanding features include simpler charging thanks to a stable stand and handy storage of the device.





TRANSPORT SAFELY

Protected storage and practical key lock



ALWAYS AT HAND

This space-saving case was specially designed to store a single NOVAFON device, as well as a power adapter/cable or your favourite attachments. Additional individual accessories can be comfortably stowed in the mesh pocket in the lid.

PRACTICAL KEY LOCK

In the workspace, on the way to a home-visit, or travelling - you can protect your device from unwanted use with the practical Travel Lock, developed exclusively for DAS NOVAFON power 2.



FOR ANIMAL THERAPISTS AND PET OWNERS

Local vibration therapy for
horses, cats and dogs

NOVAFON therapeutic treatments, which have been successfully used on humans for many decades, can also be successfully used on dogs, horses and other animals with orthopaedic or neurological disorders. Targeted local vibrations provide deep stimulation to the muscles, tissues and nerves.

Animals perceive this gentle form of therapy as very pleasant and many animal therapists report that, in addition to physiological improvements, animals also show increased psychological relaxation.

NOVAFON is proven to be effective and easy to use.

APPLICATIONS

Pain relief and regeneration from
various symptoms

HORSES

Head, neck and throat muscles
jaw joint symptoms, neck
ligament symptoms, myofascial
pain

Shoulder muscles
Myofascial pain

Back muscles
Kissing spines, sensitivity in the
saddle area, myofascial pain

Sacroiliac joint
Ileosacral joint symptoms

Hindquarters
Knee symptoms, myofascial
pain

Extremities/legs
Tendon disorders, joint
symptoms

Rideability and suppleness-
Maintain mobility and support
the warm up and cool down
phases

DOGS

Head and neck muscles
Myofascial pain

**Spine and
back muscles**

Myofascial pain,
spondylosis/spondylarthrosis,
disc prolapse, DLSS – degener-
ative lumbosacral stenosis

Extremities

Elbow dysplasia ED, hip dyspla-
sia HD, arthrosis of the joints

Quality of life Improvement of
overall well-being and relaxa-
tion

CATS

Spine and back muscles

DLSS – degenerative
lumbosacral stenosis,
spondylosis/spondylarthrosis

Muscle tension
Myofascial pain

Quality of life Improvement of
overall well-being and relaxation

Our therapists advise you personally so that you can provide optimal treatment

Our goal is for animals to be relaxed, motivated and healthy. Therefore, we support all NOVAFON users with various different services:



Application guide book and app For every NOVAFON, there are more than 100 pages to download containing basic information - in the instruction manual for horses, dogs and cats. The pictures of muscles, instructions and different areas of application support treatment. You can also find our treatment guide in the NOVAFON app - practical and fast.

Load instruction manual here:
<https://novafon.com/en/downloads>

Download app here:
<https://novafon.com/en/app>

For more than 90 years NOVAFON – ELEKTROMEDIZINISCHE GERÄTE GmbH has been manufacturing electromedical devices. With the NOVAFON brand, we are ambassadors for local vibration therapy and represent quality and reliability in this area. Our tried-and-tested therapeutic devices are made in Germany. As a manufacturer of medical devices, our quality management process is certified in accordance with DIN EN ISO 13485 and our products are governed by strict safety regulations to ensure the safety and well-being of our customers.



Certified medical device

NOVAFON is certified by TÜV SÜD as a manufacturer of medical devices. Vibration therapy has been researched for its effectiveness and the results from this research have been published in numerous studies.



Made in Germany

Our quality promise: the highest standards of workmanship and materials. Long-lasting, reliable products. Complete quality assurance.




30 DAY MONEY BACK GUARANTEE

We offer a 30-day money-back guarantee on our products. This means you can test the treatment thoroughly and risk-free.

Supported by:



on the basis of a decision
by the German Bundestag

 Products
novafon.com/en/products

NOVAFON –
Elektromedizinische Geräte GmbH
Daimlerstr. 13
71384 Weinstadt Germany

info@novafon.de | +49 7151 133970
www.novafon.com



[novafon_vet](#)